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By sample finishing we mean the features of a material surface as assessed through the micro-geometrical analysis. Roughness is the result of the unevenness (micro-geometrical errors) produced on the material surface during the manufacturing process thereof.

To characterize the roughness, the specimen profile is acquired along a given length L, called base length or sampling length. Thus, by matching the measure length and the profile detected, a line parallel to the theoretical profile can be generated so that the summary of the peaks may equal to that of the valleys. The defining parameter of a real surface is the mean roughness Ra (Sa when dealing with a 3d surface), representing the amplitude of the arithmetical mean value of the absolute values of the deviations $|y|$ of the real profile compared to the mean line. As the determination of the Ra parameter is insufficient to properly characterize a surface, it is of paramount importance to determine other values such as:

a) Skewness (Rsk) portraying the asymmetry of the density function (representing the prevalent symmetry). Rsk may be:

- 1. equal to 0 when the distribution around the mean is normal (Gaussian function);*
- 2. lesser than 0 when the distribution is negative and the profile is plain;*
- 3. greater than 0 when the distribution is positive and the profile is rich in peaks.*

b) Kurtosis (Rku), which is equivalent to the peak density of the profile (representing the indentation density). Rku may be:

- 1. equal to 3 when the distribution follows the Gaussian function;*
- 2. lesser than 3, when the curve of the function is flatter than the Gaussian function;*
- 3. greater than 3 when the curve of the function is more elongated than the Gaussian function.*

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The 2-D measures were made with the contact profilometer Form Talysurf 120, while the 3-D images and the values of Sa, Ssk and Sku were obtained recurring to the non contact pro-

DATA EMISSIONE	01/10/2014	
----------------	------------	--

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filometer Talysurf CCI 3000. The three aforementioned parameters correspond respectively to the values of Ra, Rsk and Rku determined for a surface.

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Figures 3 and 4 show the 3-D representation of the sample, while Figure 5 reproduces the 2-d profile detected with the contact profilometer. The values of the roughness parameters are reported in Tab.1.

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SLActive e SLA new batch respectively). The roughness mean value of the surface implant reported in the present report is then not superior to that of surface implants produced by Straumann.

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