	CNR IMAMOTER	
	<b>TECHNICAL REPORT</b>	Pagina 1 di 8

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*By sample finishing we mean the features of a material surface as assessed through the micro-geometrical analysis. Roughness is the result of the unevenness (micro-geometrical errors) produced on the material surface during the manufacturing process thereof.*

*To characterize the roughness, the specimen profile is acquired along a given length  $L$ , called base length or sampling length. Thus, by matching the measure length and the profile detected, a line parallel to the theoretical profile can be generated so that the summary of the peaks may equal to that of the valleys. The defining parameter of a real surface is the mean roughness  $Ra$  ( $Sa$  when dealing with a 3d surface), representing the amplitude of the arithmetical mean value of the absolute values of the deviations  $|y|$  of the real profile compared to the mean line. As the determination of the  $Ra$  parameter is insufficient to properly characterize a surface, it is of paramount importance to determine other values such as:*

*a) Skewness ( $Rsk$ ) portraying the asymmetry of the density function (representing the prevalent symmetry).  $Rsk$  may be:*

- 1. equal to 0 when the distribution around the mean is normal (Gaussian function);*
- 2. lesser than 0 when the distribution is negative and the profile is plain;*
- 3. greater than 0 when the distribution is positive and the profile is rich in peaks.*

*b) Kurtosis ( $Rku$ ), which is equivalent to the peak density of the profile (representing the indentation density).  $Rku$  may be:*

- 1. equal to 3 when the distribution follows the Gaussian function;*
- 2. lesser than 3, when the curve of the function is flatter than the Gaussian function;*
- 3. greater than 3 when the curve of the function is more elongated than the Gaussian function.*

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*The 2-D measures were made with the contact profilometer Form Talysurf 120, while the 3-D images and the values of Sa, Ssk and Sku were obtained recurring to the non contact profilometer Talysurf CCI 3000. The three aforementioned parameters correspond respectively to the values of Ra, Rsk and Rku determined for a surface.*

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

*Figures 3 and 4 show the 3-D representation of the implant, while figure 5 reproduces the 2-d profile detected with the contact profilometer. The values of the roughness parameters are reported in Tab.1.*

DATA EMISSIONE	09/01/2015	
----------------	------------	--

	<b>CNR IMAMOTER</b>	
	<b>TECHNICAL REPORT</b>	Pagina 6 di 8

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*LOT 002783 implant DM TECOM treated shows a Sa and Ra values of 1.3  $\mu\text{m}$  and 1  $\mu\text{m}$  respectively.*

*Rsk parameter is positive, while Rku value is about three; these data suggest that the surface profile exhibit a lot of peaks and few valleys.*

*Ssk is about zero and Sku about 3, indicating a quite homogeneous surface, as suggested by the low standard deviation of Sa and Ra parameters.*

*Literature data show that Ra value for double acid-etched implants ranging between 0.3 and 0.5  $\mu\text{m}$  [Ji-Hyun Kim, Yeek Herr, Jong-Hyuk Chung, Seung-II Shin, Young-Hyuk Kwon, "The effect of erbium-doped: yttrium, aluminium and garnet laser irradiation on the surface microstructure and roughness of double acid-etched implants" J. Periodontal Implant Sci. (2011) 41, 234–241], inferior to the values measured on the present implant.*

Torino, 09/01/15

Firma

DATA EMISSIONE	09/01/2015	
----------------	------------	--